

ADVANCE CARE *IP*

Breath-Activated Call Device

The Advance Care range of specialised Call Devices have been designed for those with limited movement in mind.

The Breath-Activated Call Device is used by patients with little to no motor skills and only requires a minimal amount of pressure to activate a call.

The Breath-Activated Call Device can be securely clamped to any adjacent fitting, such as a bed rail, and the flexible adjustment arm allows the device to be positioned in the most appropriate location to best suit the user's needs.

The Breath-Activated Device is UL1069 listed, allowing its use in oxygen-rich environments. The filtered mouthpiece on the unit is detachable and seamlessly replaceable.

As the Breath-Activated Device is plugged into the Nurse Call Button via the pendant socket on the front of the plate, the Call Device essentially acts as an extension of the Nurse Call Button, allowing for the easy access of a Call Point without having to physically touch the wall.

Although an extension of the Nurse Call Button, the device has its own easily identifiable alert, such as 'Breath Device Call Room 1', when it is activated.

When a call is made, a high visibility reassurance light illuminates the 'CANCEL' button on the Call Button plate and the Call Button may also 'beep' periodically (if activated), providing reassurance to the user.

The Nurse Call alert is immediately displayed at the Nurse Station and sent to associated paging devices.

The Advance Care IP Nurse Call System includes a 'Pendant Removed' alarm feature which generates an alert if the Call Device or Pendant is accidentally pulled out or incorrectly plugged in, a beneficial feature with resident safety as key.

This product integrates with most Nurse Call Systems and is a valued addition to any Advance Care IP Nurse Call installation.

Features

- Activated via breath
- Purposely designed for those with restricted movement
- Durable fully-moulded fittings
- Integrates with most Nurse Call Systems
- 6.35mm industry-standard Pendant Plug
 - Self-configured alarm notification
 - Replaceable filtered mouthpiece
 - Bed clip as standard
 - Ergonomic design
 - UL1069 listed



Code	Description
AC-BA	Breath Activated Call Device

Design Specifications

Feature	Description
Length	229mm
Weight	1250 grams
Material	Silicone
Connectors	6.35mm Phono Plug
Temperature	-20°C to 50°C
Humidity Range	0 to 90%
Rating	UL1069

Compliance

Designed to comply with:

- AS/NZS 3811 Hard Wired Patient Alarm Systems
- AS/NZS 2999 Alarm Systems for the Elderly and Other Persons at Risk
- AS/ACIF 5009 Installation Requirements for Customer Cabling
- AS/NZS 3000 Wiring Rules
- AS/NZS CISPR22 (Class B) Electromagnetic Compatibility (EMC)

© 2022 Advance Care™

All rights are reserved.

Advance Care™ reserves the right to make changes in specifications and/or discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Authorised Distributor

